

# NOT APPROVED N95 USE

Unvalved KN90  
3M 9001

"Authorized Imported, Non-NIOSH Approved  
Respirators Manufactured in China" –

[April 3rd FDA EUA Appendix A](#)

[Vertical Flatfold Earloop Video](#)

[User Instructions](#)



# NOT APPROVED N95 USE

Unvalved KN95  
3M 9502+

"Authorized Imported, Non-NIOSH Approved  
Respirators Manufactured in China"

[April 3rd FDA EUA Appendix A](#)

[Vertical Flatfold Earloop Video](#)

[User Instructions](#)



# APPROVED N95 USE

Flat Fold Masks with  
Headbands  
3M 9105

NIOSH N95 TC-84-5231

[TC-84-5231](#)

[Vflex Video](#)

[User Instructions](#)



# APPROVED N95 USE

NIOSH N95  
3M 8511

[TC-84-1299](#)

[Valved Cup FFR Video](#)

[User Instructions](#)





# APPROVED N95 USE

## Body Guard Safety Gear N95

NIOSH certified to have a filter efficiency of 95% or greater against particulate aerosols free of oil.

This mask is cleared by the FDA for use in surgical environments. Electrostatically-charged filter media provides low breathing resistance.

Latex-free, staple-free head straps for comfort and Safety-Soft foam nose cushion.

[Fastenal Product Standard](#)

